

# Hyaluronic Acid Biovivendi



# Hyaluronic Acid Biovivendi

**biovivendi®**

## What is Biovivendi®?

Biovivendi is a series of cross-linked hyaluronic acid dermal fillers that are used to treat signs of aging.

The products are developed to smooth outlines and wrinkles, shape facial contours and add volume to lips. Biovivendi is safe and effective to treat the entire face creating a more beautiful and youthful look.



# Hyaluronic Acid Biovivendi

## Hyaluronic Acid


### 1.1 What is hyaluronic acid?

Hyaluronic acid is a non-sulfated glycosaminoglycan

Distributed widely throughout the body in the connective, epithelial, and neural tissues. It is one of the chief

Components of the extracellular matrix and contributes significantly to cell proliferation and migration.

The average person has roughly 15 grams of hyaluronic acid in their body, one third of which is lost and replaced each day.



### 1.2 Hyaluronic acid is vital to:

Curb premature wrinkling

Increase skin hydration

Keep skin moist and smooth

Boost skin elasticity

Enhance collagen synthesis

Nourish healthy skin cell growth

Promote firm skin tone


Support healthier joints



# Hyaluronic Acid Biovivendi

## Biovivendi Technology

- ◇ The hyaluronic acid macromolecule forms straight lines—hyaluronan chains.
- ◇ The crosslink agent BDDE will develop hyaluronan chains into reticular formation.
- ◇ Using advanced technology, reticular formation molecule become three-dimensional matrix, which will improve both viscosity and elasticity of the gel.
- ◇ After injected, the HA filler locks water and expands instantly. Then the soft tissue depression is filled.



# Hyaluronic Acid Biovivendi



## Biovivendi Finelines

Finelines is mainly used for wrinkles around eyes, glabellar wrinkles, crow's feet lacrimal sulcus and nasolabial etc; usually be injected into shallow dermis.

### SPECIAL FEATURES

- ◆ Non-animal origin: sodium hyaluronate synthesized by bio-fermentation
- ◆ Cross-linking: with butanediol diglycidyl ether (BDDE), the cross-linking agent universally recognized as the least toxic of all, which ensures that the risk of allergies is minimal.
- ◆ Dual-action: high-concentration filler: 20mg/ml of cross-linked hyaluronic acid with bio-stimulating free hyaluronic acid



## Biovivendi Finelines

Content	1ml, 2ml
Concentration	20mg/ml
Reticulating Agent	BDDE
Particle Size	0.10-0.15mm

# Hyaluronic Acid Biovivendi

## Biovivendi DERM

Derm is mainly used for lip wrinkles, brow wrinkles, correct of aging atrophy, dripping lips, sag cheek, facial asymmetry and fill earlobe; usually be injected into middle or deep dermis.

### SPECIAL FEATURES

- ◆ Non-animal origin: sodium hyaluronate synthesized by bio-fermentation
- ◆ Cross-linking: with butanediol diglycidyl ether (BDDE), the cross-linking agent universally recognized as the least toxic of all, which ensures that the risk of allergies is minimal.
- ◆ Dual-action: high-concentration filler: 24mg/ml of cross-linked hyaluronic acid with bio-stimulating free hyaluronic acid



### Biovivendi DERM

Content	1ml, 2ml
Concentration	24mg/ml
Reticulating Agent	BDDE
Particle Size	0.15-0.28mm



# Hyaluronic Acid Biovivendi

## Biovivendi Deep

Deep is mainly used for filling or shaping of face or nose, plumping the lips, forehead, chin etc. Nalsal plastic surgery effect: improve the poor nasal shape, low nose, saddle nose and nasal tip. Usually be injected in deep dermis.

### SPECIAL FEATURES

- ◇ Non-animal origin: sodium hyaluronate synthesized by bio-fermentation
- ◇ Cross-linking: with butanediol diglycidyl ether (BDDE), the cross-linking agent universally recognized as the least toxic of all, which ensures that the risk of allergies is minimal.
- ◇ Dual-action: high-concentration filler: 24mg/ml of cross-linked hyaluronic acid with bio-stimulating free hyaluronic acid



### Biovivendi Deep

Content	1ml, 2ml
Concentration	24mg/ml
Reticulating Agent	BDDE
Particle Size	0.28-0.50mm


# Hyaluronic Acid Biovivendi

## Advantages and Distinguishing Features:

- ◆ The treatment is immediate, since it requires only a few minutes and can be carried out straight after a simple medical consultation
- ◆ The application is virtually painless and generally requires no anesthesia of any sort.
- ◆ The glass syringe is very easy to handle, practical for application, allow considerable control over the amount injected and give an excellent finished filling effect.
- ◆ The result of treatment are long-lasting but not permanent; indeed, over time natural reabsorption of the pure hyaluronic acid occurs.

### Main Treatment Areas

- Brow furrows
- Cheek augmentation
- Lips contour
- Lips augmentation
- Redefinition of the facial contour
- Temples
- Periorbital lines
- Nose augmentation
- Nasolabial folds
- Chin augmentation





# Hyaluronic Acid Biovivendi



Good Effectiveness



Before and After



# Hyaluronic Acid Biovivendi

Good Effectiveness



Before and After



# Our Company

Biovivendi production and trading company was established in April 2011. This hyaluronic acid syringes come from a high and new technology enterprise integrating the R&D, production and sales of medical biomaterial. With well-qualified staff, thorough quality management system and first-class facilities, the company ensures the continuous production of injection products that confirm to the laws and regulations of the State Food and Drug Administration and EU MDD.

